

# ADVANCED FALL&VITALS ALERT SYSTEM

## PROJECT REPORT

MALAVIKA KR

25BOE10137

# INTRODUCTION

This project outlines the development of a wearable, real-time patient monitoring system designed to automatically detect and respond to medical emergencies, specifically focusing on patient falls.

The system utilizes an accelerometer to constantly analyze the patient's motion and a heart rate sensor to assess their physiological state. The core objective is to minimize the time between a critical event (a fall) and the commencement of emergency aid.

# PROBLEM STATEMENT

The core problem is the critical delay and poor reliability inherent in current fall alert systems for at-risk individuals.

Falls are major injury causes, but existing systems either fail to notify if the patient is incapacitated (the long-lie problem) or generate too many false alarms, wasting caregiver time.

# FUNCTIONAL REQUIREMENTS

Detect Impact: Trigger on  $\geq 2.5G$

Confirm Fall:  $\geq 10$  seconds stillness ( $\leq 0.2G$ ).

Check Pulse: Wait 3s. Alert if  $< 50$  or  $> 20$  BPM.

Notify: Send alert on confirmation/  
critical pulse

set **is\_patient\_safe** to **FALSE**.

# NON FUNCTIONAL REQUIREMENTS

Performance: Real-time data processing ( $< 100$  ms latency).

Reliability: High accuracy ( $\geq 95\%$  TPR,  $\leq 1\%$  FPR).

Availability: High uptime (99.9%).

Security: Data encryption (HIPAA/GDPR) and access authentication.

Constraint: Battery life  $\geq 24$  hours.

# SYSTEM ARCHITECTURE


Wearable: Sensors (Accel + Pulse) to  
MCU (Logic  $\geq$  2.5G / 10sStillness).

Gateway: Relays Wearable data (BLE) to Cloud (Wi-Fi/Cellular).

Cloud Server: Validates Alerts, Logs Data, Sends Notifications.

UI: Caregiver Dashboard for Alerts and Acknowledgment.

# DESIGN DIAGRAMS



# IMPLEMENTATION DETAILS

Fall Confirmation (Accelerometer): A high impact followed by 10 seconds of stillness confirms the fall and triggers a primary alert.

Condition Check (Pulse): 3seconds after the fall, the patient's heart rate is checked. An alert is critically escalated if the pulse is outside the 50 to 120 BPM



```
❏ # --- Constants ---  
FALL_G = 2.5  
STILL_G = 0.2  
PULSE_LOW = 50  
PULSE_HIGH = 120  
  
# --- Simple Fall + Pulse Check ---  
g = float(input("Enter G-force: "))  
  
if g > FALL_G:  
    print(" Impact detected!")  
  
    still = float(input("Enter stillness G-force: "))  
    if still < STILL_G:  
        print("Fall Confirmed!")  
  
        bpm = int(input("Enter Heart Rate: "))  
        if bpm < PULSE_LOW or bpm > PULSE_HIGH:  
            print("Critical Pulse Alert!")  
        else:  
            print("Pulse Normal")  
    else:  
        print(" No Fall. Patient moved.")  
  
else:  
    print("Normal Activity")
```

# TESTING APPROACH

Unit Testing: Verify individual constants work

Integration Testing: Test the full logic flow using simulated data.

System Testing: Test real-world reliability using physical scenarios (deliberate falls, ADL false alarms) and check alert latency

# CHALLENGES FACED

- False Alarms: High rate of False Positives (e.g., sitting heavily) and False Negatives (missing slow falls).
- Calibration: Fixed thresholds (2.5G) don't adapt to individual patients or environments.
- Power: Continuous monitoring of G-force and pulse is a major battery drain.
- Data Quality: Motion artifacts corrupt heart rate data after a fall.
- Compliance: Device must be comfortable for 24/7 wear to be effective.

# LEARNING AND TAKEAWAY

- Learning: Simple, static thresholds (like 2.5G for impact) are efficient but fundamentally limited. They cannot distinguish between a real, dangerous fall and vigorous everyday activities.
- Takeaway: Robust safety systems require adaptive algorithms or machine learning to minimize false alarms (reducing caregiver "alarm fatigue") while ensuring high sensitivity to genuine falls.

# FUTURE ENHANCEMENTS

- AI/ML: Implement Adaptive Thresholds to reduce false alarms.
- Sensor Fusion: Add Gyroscope and Ambient Sensors for better fall confirmation.
- Optimization: Use Edge Computing for better battery life and integrate with telemedicine.

# REFERENCES

1. Academic Papers/Articles Referred
2. World Health Organization. "Patient safety." WHO Fact Sheet.
3. Oliver D, Healey F, Haines TP. "Hospital Fall Prevention: A Systematic Review of Implementation, Components, Adherence, and Effectiveness." J Am Geriatr Soc. 2013-03-24.  
<https://pmc.ncbi.nlm.nih.gov/articles/PMC3670303/>